

Company Profile

Shandong HonHai Industrial Co. Ltd.

We started to produce neutral pharmaceutical glass tube from the year 1998, and become the leading manufacture in China of this important raw material for tubular glass containers. In the following years, we rapidly expanded our business to nearly cover all pharmaceutical glass package. Now we make or supply the follow items.

- Neutral pharmaceutical glass tube UPS TYPE I
(COE 5.0 & COE 7.0)
- Tubular glass vial and ampoule
- Moulded injection glass vial
- Moulded infusion glass bottle
- USP TYPE I glass
- Cartridge & prefillable syringe
- Blood collect tube
- Glass dropper
- Butyl rubber stopper
- Aluminum plastic cap










This section will focus on one important branch: Amber glass vial and bottle.



Ref.	ml	ml	g			mm
DZK-8	20	26	29	ISO/SFDA	58	32
DZK-27	30	38	35	ISO/SFDA	62.8	36
DZK-28	50	60	50	ISO/SFDA	73	41
DZK-9	100	119	89	ISO/SFDA	94.5	51.6
DZK-114	100	117	92	ISO/SFDA	94.5	51
DZK-115	50	58	52	ISO/SFDA	73	42
DZK-125	15	17	22	ISO/SFDA	58.8	26.5










						
DZK-13	10	12.5	27.5	DIN18	58	24.8
DZK-126	15	20	35	DIN18	65	28.5
DZK-15	20	22.5	35	DIN18	71.5	28.7
DZK-3	30	35	48	DIN18	78.8	33
DZK-16	50	57	55	DIN18	92.2	37.2










Ref.	ml	ml	g	mm	mm	mm
DZK-63	5	8	23	20mm	55	22
DZK-64	10	13.5	26	20mm	65	24
DZK-65	15	18.5	30	20mm	74	25.5
DZK-66	20	23.5	32	20mm	80	26.5
DZK-67	30	35	50	20mm	96	30
DZK-68	50	55.5	56	20mm	102	35
DZK-69	100	105	95	20mm	124	43










						
DZK-53	14.2	16.4	25.5	18mm-400	68.26	25
DZK-54	31.3	36.2	44	20mm-400	78.58	32.8
DZK-55	60.8	63.8	58	20mm-400	93.66	38.6



						
DZK-84	15	18.5	32	PP18	67.5	27
DZK-86	30	34.5	45	20	78.7	33
DZK-57	30	35	50	PP22	84.5	32.5
DZK-61	60	65	74.8	20.15	93.85	40.4
DZK-83	60	66	71	24.6	96	39.7
DZK-70	60	71	68	PP28	97.3	41
DZK-71	100	110	95	PP28	101.6	49



						
DZK-39	20	26.5	45	PP22	80.5	28.5
DZK-38	20	27	47	PP22	81.4	29
DZK-6	30	36	50	PP25	88.5	30
DZK-7	30	35	50	PP22	89.5	30
DZK-4	50	58	82	PP28	107.2	35










Ref.	ml	ml	g		mm	mm
DZK-72	100	115.5	109	27.5	115	46
DZK-1	100	110	109	PP 28	120.7	44.5
DZK-58	100	113	106	PP-28	120.7	44.3
DZK-59	100	110	113	GRP-0106-28	120	44.6
DZK-32		138	128	27.4	140	46










Ref.	ml	ml	g			mm
DZK-29	100	127	114	27.6	122	47
DZK-40	100	119	117	27.8	120.5	44.5
DZK-52	100	109.5	120	27.6	121.5	44.3
DZK-30		136	115	27.9	130	48.8
DZK-60	153	160	145	27.7	136	52











						
DZK-75	10	12.5	26	12.87	57	32
DZK-76	20	23	45	17.5	69.1	38.5
DZK-77	30	33	50	17.5	75.5	42
DZK-78	50	53	74	17.5	86	49
DZK-79	100	108	115	17.5	106	60










						
DZK-34		10	22	17.7	50.5	32
DZK-35		10	23.5	18	51.5	30
DZK-51	6	11.5	24.5	18	50	30.5
DZK-62		10.5	24	18.3	51	28.5










						
DZK-92		16	28	17.4	67.3	25
DZK-93		33.5	38	19.6	78	30.4
DZK-94		66	65	19.6	94	40
DZK-95	120	126	108	21.5	113.3	48.8



						
DZK-80	60	67.5	60	27.7	94.5	38.6
DZK-81	125	138	103	27.9	94.5	38.6
DZK-88	100	108	72	27.7	104.8	46.1
DZK-87	30	35.5	40	27.7	67	35
DZK-96	10	12	26	18.1	54.8	25
DZK-103	15	19	30	21.9	63	27.5










						
DZK-106	10	15	36.5	17.5	58.5	33.3
DZK-107	20	24	51.5	17.5	65.7	39
DZK-108	30	34	60	17.5	72.7	41
DZK-109	50	56	75	17.5	90	45.5
DZK-110	100	108	115	17.5	113.5	55



Ref.	ml	ml	g		mm	mm
DZK-104	100	110	108	28	115.4	47
DZK-122	100	115	110	28	125.25	42.9
DZK-117	100	105	144	25	137.6	53
DZK-116	100	110	125	28	138	43.2
DZK-111	100	110	115	28	138	42.4
DZK-11	100	110	118	28	126	44.6
DZK-99		61	65	28	94	36.5



						
DZK-97	60	65	85	27	80.5	44
DZK-98	30	36	58	27	68	36.7
DZK-118	60	69	73	33	75.5	44.5
DZK-100	100	114	92	38.5	87.35	50
DZK-113	40	50	51	28	79.5	36.5
DZK-33		9.5	22	21.5	43	25.6
DZK-37		33	65	28.8	66.2	43.5